



## **EURO PLUS Infection Hip Spacer**



**Your partner for implants and  
trauma products**

## **EURO PLUS Infection Hip Spacer**

In total hip arthroplasty, infected periprosthetic femoral fractures are some of the most complex and significant complications.

By using EURO PLUS Infection Hip Spacer, surgeons are able to create intraoperative, temporary cement spacers corresponding to modular cemented stem and femoral head for patients undergoing a 2- Steps revision due to infection.

The clearly shaped and smooth surface allows the patient a certain mobility during the implant period. Temporary antibioticloaded cement spacer facilitates the extermination of infections and fracture stabilization.


The EURO PLUS Infection Hip Spacer Molds Instrument can be sterilized for reuse.

We recommend to use bone cement with standard viscosity to making Hip cement spacers.



## EURO PLUS Infection Hip Spacer

EURO PLUS Infection Hip Spacer - Stem		
Ref. Number	Stem Size	Stem Length L (mm)
EUTS01P001N	0	127
EUTS01P002N	1	131
EUST01P003N	2	135
EUST01P004N	3	139
EUST01P005N	4	143
EUST01P006N	5	147
EUST01P007N	6	151
EUST01P008N	7	155
EUST01P009N	8	159
EUST01P010N	9	163



EURO PLUS Infection Hip Spacer - Head	
Ref. Number	Head Size (Diam. mm)
EUTS02P046N	46
EUST02P048N	48
EUST02P050N	50
EUST02P052N	52
EUST02P054N	54
EUST02P056N	56
EUST02P058N	58
EUST02P060N	60


EURO PLUS Infection Hip Spacer - Foot Stand	
Ref. Number	
EUTS03P001N	




EURO PLUS Infection Hip Spacer - Keel	
Ref. Number:	EUTS04P000N- EUTS04P009N
10 sizes, it depends on the the length of stem	

**EURO PLUS Infection Hip Spacer  
Instrument Set**  
**EK0900A001N**

**Tray 1**



**Tray 2**



**Tray 3**  
**( Optional: Bone Cement Removal Instruments )**







**EURO PLUS Infection Hip Spacer  
Instrument Set ( Tray 1)**

**EK0900A001N**

No	Ref. Number	Description
1	EUTS02P046N	EURO PLUS Infection Spacer Head - 46
2	EUTS02P048N	EURO PLUS Infection Spacer Head - 48
3	EUTS02P050N	EURO PLUS Infection Spacer Head - 50
4	EUTS02P052N	EURO PLUS Infection Spacer Head - 52
5	EUTS02P054N	EURO PLUS Infection Spacer Head - 54
6	EUTS02P056N	EURO PLUS Infection Spacer Head - 56
7	EUTS02P058N	EURO PLUS Infection Spacer Head - 58
8	EUTS02P060N	EURO PLUS Infection Spacer Head - 60
9	EUTS01P001N	EURO PLUS Infection Spacer Stem - Size 0
10	EUTS01P002N	EURO PLUS Infection Spacer Stem - Size 1
11	EUTS01P003N	EURO PLUS Infection Spacer Stem - Size 2
12	EUTS01P004N	EURO PLUS Infection Spacer Stem - Size 3
13	EUTS01P005N	EURO PLUS Infection Spacer Stem - Size 4
14	EUTS01P006N	EURO PLUS Infection Spacer Stem - Size 5
15	EUTS01P007N	EURO PLUS Infection Spacer Stem - Size 6
16	EUTS01P008N	EURO PLUS Infection Spacer Stem - Size 7
17	EUTS01P009N	EURO PLUS Infection Spacer Stem - Size 8
18	EUTS01P010N	EURO PLUS Infection Spacer Stem - Size 9
19	EUTS03P001N	EURO PLUS FOOT STAND



**EURO PLUS Infection Hip Spacer  
Instrument Set ( Tray 2 )**

**EK0900A001N**

No	Ref. Number	Description
1	EE9937N046N	CEMENT SPACER MOLD FOR UNIPOLAR HEAD / Ø46
2	EE9937N048N	CEMENT SPACER MOLD FOR UNIPOLAR HEAD / Ø48
3	EE9937N050N	CEMENT SPACER MOLD FOR UNIPOLAR HEAD / Ø50
4	EE9937N052N	CEMENT SPACER MOLD FOR UNIPOLAR HEAD / Ø52
5	EE9937N054N	CEMENT SPACER MOLD FOR UNIPOLAR HEAD / Ø54
6	EE9937N056N	CEMENT SPACER MOLD FOR UNIPOLAR HEAD / Ø56
7	EE9937N058N	CEMENT SPACER MOLD FOR UNIPOLAR HEAD / Ø58
8	EE9937N060N	CEMENT SPACER MOLD FOR UNIPOLAR HEAD / Ø60
9	EE1039S000N	EUROFIT Primary Cementless Stem Rasp / 0
10	EE1039S001N	EUROFIT Primary Cementless Stem Rasp / 1
11	EE1039S002N	EUROFIT Primary Cementless Stem Rasp / 2
12	EE1039S003N	EUROFIT Primary Cementless Stem Rasp / 3
13	EE1039S004N	EUROFIT Primary Cementless Stem Rasp / 4
14	EE1039S005N	EUROFIT Primary Cementless Stem Rasp / 5
15	EE1039S006N	EUROFIT Primary Cementless Stem Rasp / 6
16	EE1039S007N	EUROFIT Primary Cementless Stem Rasp / 7
17	EE1039S008N	EUROFIT Primary Cementless Stem Rasp / 8
18	EE1039S009N	EUROFIT Primary Cementless Stem Rasp / 9

**EURO PLUS Infection Hip Spacer  
Instrument Set ( Tray 3 )**

**EK0900A001N**

No	Ref. Number	Description
1	EUB084R001N	Curette, 10mm
2	EUB093R002N	Osteotome, straight 9mm, 400mm
3	EUB087R003N	Pituitary 5,3x20mm, 260mm
4	EUB089R004N	Tap 14mm
5	EUB080R005N	Tap 10mm
6	EUB081R006N	Tap 12mm
7	EUB083R007N	Curette 8mm
8	EUB113R008N	Slide Hammer 680g
9	EUB091R009N	Osteotome, curved 9mm, 400mm
10	EUB066R010N	Osteotome, half edge 12mm, 400mm



## HOW TO CONTACT US

Euromed Implants GmbH  
Wilhadi-Kirchhof 13  
21682 Stade, Germany

Phone: +49 (0) 4141 7897 290  
Fax: +49 (0) 4141 7782 700

[office@euromed-implants.de](mailto:office@euromed-implants.de)  
[inquiry@euromed-implants.de](mailto:inquiry@euromed-implants.de)  
[www.euromed-implants.de](http://www.euromed-implants.de)